IKA® lab disc

<i>BETRIEBSANLEITUNG</i>	<i>D</i>	<i>3</i>
<i>OPERATING INSTRUCTIONS</i>	<i>GB</i>	<i>6</i>
<i>MODE D'EMPLOI</i>	<i>F</i>	<i>9</i>
<i>INSTRUCCIONES DE MAEJO</i>	<i>E</i>	<i>12</i>
<i>HANDLEIDING</i>	<i>NL</i>	<i>15</i>
<i>ISTRUZIONI PER L'USO</i>	<i>I</i>	<i>18</i>
<i>DRIFTSANVISNING</i>	<i>S</i>	<i>21</i>
<i>DRIFTSINSTRUKS</i>	<i>DK</i>	<i>24</i>
<i>DRIFTSVEJLEDNING</i>	<i>N</i>	<i>27</i>
<i>KÄYTTÖOHJE</i>	<i>SF</i>	<i>30</i>
<i>INSTRUÇÕES DE SERVIÇO</i>	<i>P</i>	<i>36</i>

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D

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Ex -Schutz:  II 3G EEx nC IIB T6

Contens

	Page
Safety instructions	6
Correct use	6
Unpacking	6
Commissioning	6
Maintenance and cleaning	7
Guarantee	7
Technical data	8

Safety instructions

This instrument is designed for use in Explosion Hazard Zone 2 areas. Do not use in Explosion Hazard Zone 1 areas.

If you are working in explosion hazard environments, do not use plastic beakers. Plastic beakers can create static electrical charges and can ignite explosive substances. The power supply plug must be connected up outside the explosion hazard zone. The housing must be earthed. To earth the instrument, use a earthing cable with a commercial 6.3 quick-connect terminal and plug this into the socket on the rear of the instrument (see figure).

Damage to the protective membrane sheet will reduce the IP and explosion hazard zone rating of the instrument.

Hazard warning – Magnetism. The instrument creates a magnetic field. Items which are sensitive to magnetism (e.g. data storage devices, pacemakers...) must be kept away from the instrument.

Correct use

The **IKA® lab disc** is a magnetic stirrer designed for mixing up to 800 ml of thinly fluid media with magnetic bars up to 30 mm long. Stirring speed can be adjusted in the range from 15 to 1500 rpm.

Unpacking

Please unpack the device carefully and check it for damage. It is important that any transport damage is detected when the device is unpacked. If necessary, an inventory of the damage should be made immediately (postal service, railway, haulage company).

As supplied, the package contains:

The **IKA® lab disc**, a plug-in power supply unit and an instruction manual.

Commissioning

Plug the instrument into a power supply socket. (For the power supply specifications, see the section "Technical Data").

Place the magnet bar and the medium to be mixed in a container which is transparent to magnetic fields. Set the stirrer speed to minimum. To do this, switch on the instrument by pressing the  button and then press the  button until the flashing LED between the minus and plus buttons stops flashing and lights permanently. You can now place the container on the mark in the centre of the set-up surface.

Now press the \ominus button and the \oplus button to set the speed you require. The LED between the buttons will flash during setting.

To select the reverse mode and invert the direction of stirring at regular intervals, press the \odot button. In the reverse mode the stirring direction will be inverted every 30 seconds; there will be a pause of 5 seconds between each changeover. The LED above the \odot button will flash when the reverse mode is active.

The speed setting and reverse mode selection will remain unchanged when the instrument is unplugged.

Caution: At some speeds, irregularities in the shape of the stirring bar and of the bottom of the container can cause the bar to move irregularly during stirring. Increase or decrease the speed until the stirring bar begins to rotate evenly.

Maintenance and cleaning

The **IKA® lab disc** is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate. Use only cleansing agents which have been approved by IKA to clean IKA devices.

To remove use:

Dyes	isopropyl alcohol
Constructions materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning. Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Technical data

Plug-in power supply unit:

Rated voltage: **V AC** 100 - 240

Rated frequency: **Hz** 50 - 60

Output voltage: **VDC** 12

Absorbed power: **W** 5

Output power: **W** min. 3

Magnetic stirrer: coil drive 4-pole

Speed range: **1/min** 15 - 1500

Stirred quantity (max.): **ml** 500 (at 1500 rpm with
25mm stirring bar)

Operating time: **%** 100

Protection class: II

Pollution degree: 2

Overvoltage category: II

EN 60529 protection class: IP 65

Ambient humidity (relative) **%** 80

Ambient temperature: **°C** 5 to 40 (in the operation)

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