

## Data Sheet

Mixing / Magnetic Stirrers



### topolino mobil

Extremely convenient magnetic mini-stirrer for mixing quantities up to 250 ml.

- Durable, brushless motor
- Continuously adjustable speed range
- High magnetic adhesion
- Portable operation with long operating time (8-12 h)
- Short charging time (2-3 h)
- Standard replaceable AA rechargeable batteries
- Optional power mode:
  - a) Mains-free with standard batteries
  - b) With supplied mains adapter (without batteries)
  - c) Combined mains/battery operation (with batteries fitted)

Accessories: TRIKA® 25 Magnetic stirring bar, RSE Stirring bar remover, IKAFLON® 10 Magnetic stirring bar, IKAFLON® 30 Magnetic stirring bar, IKAFLON® 20 Magnetic stirring bar, IKAFLON® 25 Magnetic stirring bar, IKAFLON® 15 Magnetic stirring bar

#### Technical Data

Number of stirring positions	1
Speed range [1/min]	300 - 1800
Set-up plate material	synthetic (PP)
Protection class according to DIN EN 60529	IP 20
Voltage [V]	100 - 240
Frequency [Hz]	50/60
DC Voltage [V=]	9
Power input [mA]	1000
<b>Ident. No.</b>	<b>3381301</b>