

## IKA®-PET 50µl 1ch



Light, high quality air cushion flask lift pipette for science, research and routine work in the field of liquid handling. The IKA® LAB Pet Pipette, which was developed in accordance with the latest ergonomic findings, includes the following features:

- Pipette with certificate of conformity
- Maximum accuracy and precision of pipette and tip system
- High quality casing and flask material
- Minimum manpower required even for blow out and tip ejection
- Easy to operate
- Easy to maintain thanks to the assembly tool provided
- Adjustable even for liquids with a different density to water
- Two year warranty

Technical Data			
Design [-]	piston stroke pipette - air cushion	Handling [-]	mechanical
Volume [-]	fix	Volume display [-]	no
Number of channels [-]	1	Nominal volume [µL]	50
Accuracy: Systematic error es [±%]	0,7	Precision: Random error CV [≤%]	0,3
Certificat of conformity according to DIN 8655 [-]	yes	Readjustable oder can be readjusted [-]	yes
Blow-Out Function [-]	yes	Tip ejection [-]	yes
can be autoclaved [-]	yes	Dimensions (W x H x D) [mm]	28 x 206 x 52
Weight [kg]	0,075		