

MF 10 basic Microfine grinder drive



Continuously operating universal grinder.

- Powerful drive
- Easy to clean working surface made of stainless steel
- Two different grinding heads can be attached to the drive
- Heads are easily changeable
- Grinding heads not incl. with delivery

Technical Data			
Process type [-]	inline	Operating principle [-]	cutting/impact
Motor rating input [W]	1000	Motor rating output [W]	500
Speed range [1/min]	3000 - 6500	Power-on time ON [min]	120
Power-on time OFF [min]	30	Milling chamber, can be cooled with water [-]	no
Mill feed can be cooled in milling chamber with dry ice [-]	yes	Mill feed can be cooled in milling chamber with liquid nitrogen [-]	no
Dimensions (W x H x D) [mm]	320 x 380 x 300	Weight [kg]	9.7
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 22	RS 232 interface [-]	no
Analog output [-]	no	Voltage [V]	230 / 115
Frequency [Hz]	50/60	Power input [W]	1000