
IKA Service-Manual

Disperser / Dispergiergerät T 50 basic ULTRA-TURRAX

No. 29531.00



IKA®

IKA – Werke GmbH & Co. KG
Janke & Kunkel-Str. 10
D - 79219 Staufen

Phone: 0049 (0)7633 831 139 / 202
Fax: 0049 (0) 7633 500 212
E-Mail: service@ika.de
Internet: www.ika.de

Type: **T 50 basic ULTRA TURRAX**Item No.:**2953100**

Voltage: 230V~ 50/60Hz

IKA® Service-Manual

Date: 25-APRIL-2007

SZ: 2953100SZ_a

Pos	Item No	Bezeichnung	Description	Technical data	Bemerkung	Remark
1	2952800	UNIVERSALMOTOR	UNIVERSAL MOTOR	n10000V 230W1100,0ac		
2	2952100	FLANSCH	FLANGE			
3	2952600	FRONTSCHILD/FOLIE	FRONT FOIL	66,0x 15,0		
4	2952700	DREHZAHLTABELLE/FOLIE	SPEED CHART	24,0x 40,0		
5	1637401	KLEMMSCHRAUBE	CLAMPING SCREW			
6	2952400	KUPPLUNG	COUPLING			
7	2600800	AUSLEGER	BRACKET	D 16xL 220		
8	0864600	ZYL.SCHRAUBE	SOCKET HEAD CAP SCREW	M 5,0x 20,0 8.8-A		
10	2746600	WARNSCHILD	DANGER SIGN	SL 12,5		
3001	3032300	ZULEITUNGSKABEL	FEED CABLE			
3002	3032400	KABELVERSCHRAUBUNG	CABLE CONNECTION			
3003	3032500	ELEKTRONIK	ELECTRONIK PCB			
3004	3032600	KOHLE	MOTOR CARBON			
3005	3032700	RILLENKUGELLAGER	BALL BEARING	10x 30x 9 6200-2RSC3		
3006	3032800	RILLENKUGELLAGER	BALL BEARING	12x 32x10 6201-2RSC3		
3007	3032900	LAGERGUMMI	STOCK RUBBER			
3008	3033000	RINGMAGNET	RING MAGNET			
3009	3033100	AUFNAHMEBOLZEN	RECORDING BOLT			
3010	3033200	DISTANZHÜLSE	DISTANCE SLEEVE			
3011	3033300	ANKER MIT LUEFTER	ANCHOR			
3012	3033400	KOHLLENHALTER	CARBON HOLDER			
3013	3033500	SCHALTER	SWITCH			
3014	3033600	KONDENSATOR	CAPACITOR			
3015	3240700	MOTORGEHÄUSE	MOTOR HOUSING			
3016	3247900	SCHALTERGEHÄUSE	SWITCH HOUSING			
3017	3248000	GITTER	TRELLIS			

