

COLOR SQUID IKAMAG® Curly



Small magnetic stirrers without heating, in a variety of motifs.

- Recyclable materials
- Very good chemical resistance due to glass top and synthetic bottom made of Hytrel®
- Electronically controlled motor with infinitely variable speed from 0 - 1,500 rpm

| Technical Data | | | |
|------------------------------------------------|----------------|---------------------------------------------------------------------|-----------|
| Number of stirring positions [-] | 1 | Stirring quantity max. per stirring position (H ₂ O) [l] | 0.8 |
| Motor rating input [W] | 2 | Motor rating output [W] | 1 |
| Speed range [1/min] | 0 - 1500 | Stirring bar length max. [mm] | 30 |
| Set-up plate material [-] | glass | Set-up plate dimensions [mm] | Ø 110 |
| Dimensions (W x H x D) [mm] | 130 x 50 x 135 | Weight [kg] | 0.48 |
| Permissible ambient temperature [°C] | 5 - 40 | Permissible relative moisture [%] | 80 |
| Protection class according to DIN EN 60529 [-] | IP 54 | Voltage [V] | 100 - 240 |
| Frequency [Hz] | 50/60 | Power input [W] | 4 |