

Product Data Sheet



COLOR SQUID IKAMAG® Harry's note



Small magnetic stirrers without heating, in a variety of motifs.

- Recyclable materials
- Very good chemical resistance due to glass top and synthetic bottom made of HytreI®
- Electronically controlled motor with infinitely variable speed from 0 - 1,500 rpm

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H2O) [l]	0.8
Motor rating input [W]	2	Motor rating output [W]	1
Speed range [1/min]	0 - 1500	Stirring bar length max. [mm]	30
Set-up plate material [-]	glass	Set-up plate dimensions [mm]	Ø 110
Dimensions (W x H x D) [mm]	130 x 50 x 135	Weight [kg]	0.48
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 54	Voltage [V]	100 - 240
Frequency [Hz]	50/60	Power input [W]	4