

## KMO 2 basic IKAMAG®



Small, powerful magnetic stirrer without heating.

- Strong magnetic field
- Motor with optoelectronic speed control
- Infinitely variable speed from 0 - 1,100 rpm
- Stainless steel casing facilitates cleaning and sterilization
- Includes M 10 thread for H 16 V support rod

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	5
Stirring quantity max. (H <sub>2</sub> O) [l]	5	Motor rating input [W]	14
Motor rating output [W]	4	Speed range [1/min]	0 - 1100
Stirring bar length max. [mm]	50	Set-up plate material [-]	stainless steel 1.4301
Set-up plate dimensions [mm]	140 x 120	Dimensions (W x H x D) [mm]	140 x 75 x 200
Weight [kg]	1.4	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 21
Voltage [V]	230 / 115 / 100	Frequency [Hz]	50/60
Power input [W]	20		