

Midi MR 1 digital IKAMAG®



- Powerful magnetic stirrer without heating.
- Flat, sturdy stainless steel casing
 - Non-locking motor
 - Infinitely variable speed
 - Digital LED speed display
 - Timer (0 - 56 min) or continuous operation
 - For stirring quantities up to 50 liters (H₂O)

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H ₂ O) [l]	50
Motor rating input [W]	70	Motor rating output [W]	19
Speed display [-]	LCD	Speed range [1/min]	0 - 1000
Stirring bar length max. [mm]	80	Set-up plate material [-]	stainless steel 1.4301
Set-up plate dimensions [mm]	350 x 350	Dimensions (W x H x D) [mm]	360 x 110 x 430
Weight [kg]	9	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 21
Voltage [V]	230 / 100 - 120	Frequency [Hz]	50/60
Power input [W]	70		