

REO basic C IKAMAG®



Classic magnetic stirrer without heating new designed. Outstanding chemical resistance due to the white coated heating plate.

- Non-locking, electronically controlled motor
- Constant speed even during changes in load
- Infinitely variable speed

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H ₂ O) [l]	20
Stirring quantity max. (H ₂ O) [l]	20	Motor rating input [W]	12
Motor rating output [W]	5	Speed range [1/min]	0 - 1500
Stirring bar length max. [mm]	80	Set-up plate material [-]	technical enamel
Set-up plate dimensions [mm]	Ø 135	Dimensions (W x H x D) [mm]	160 x 97 x 280
Weight [kg]	2.4	Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80	Protection class according to DIN EN 60529 [-]	IP 42
Voltage [V]	220 - 230 / 115 / 100	Frequency [Hz]	50/60
Power input [W]	20		