

## RO 10 power IKAMAG®



Multi-position magnetic stirrer with 10 stirring positions, without heating. The stainless steel surface covers the unit allowing easy cleaning and providing protection against the penetration of liquids.

- Optimum use of laboratory space
- Including removable PUR cover

Technical Data			
Number of stirring positions [-]	10	Stirring quantity max. per stirring position (H <sub>2</sub> O) [l]	0.4
Stirring position distance [mm]	90	Speed deviation of individual stirring positions [%]	5
Stirring quantity max. (H <sub>2</sub> O) [l]	4	Motor rating input [W]	14.4
Motor rating output [W]	3.6	Speed range [1/min]	0 - 1100
Stirring bar length max. [mm]	30	Speed control [-]	scale 1-10
Set-up plate material [-]	stainless steel 1.4301	Set-up plate dimensions [mm]	180 x 450
Dimensions (W x H x D) [mm]	182 x 65 x 552	Weight [kg]	3.2
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 42	Voltage [V]	230 / 115
Frequency [Hz]	50/60	Power input [W]	20