

lab disc (SUNNY SIDE UP)



Ultra-flat compact magnetic stirrer, guaranteed with modern magnet coil technology. Wear-free drive with no moving parts. The lab disc can reverse direction of rotation automatically every 30 seconds to ensure better mixing.

- Explosion hazard
- Rotation direction is reversible
- High IP protection class (IP 65)
- Wear-free
- Set-up plate and casing made from chemically resistant materials
- Slip-proof, safe stand

Technical Data			
Number of stirring positions [-]	1	Stirring quantity max. per stirring position (H ₂ O) [l]	0.8
Motor rating input [W]	5	Motor rating output [W]	3
Speed range [1/min]	15 - 1500	Stirring bar length max. [mm]	25
Set-up plate material [-]	polyester	Set-up plate dimensions [mm]	Ø 90
Dimensions (W x H x D) [mm]	114 x 12 x 161	Weight [kg]	0.3
Permissible ambient temperature [°C]	5 - 40	Permissible relative moisture [%]	80
Protection class according to DIN EN 60529 [-]	IP 65	Voltage [V]	100 - 240
Frequency [Hz]	50/60	Power input [W]	5