

HS 260 basic



- Compact, flat shaker with ideal swivel motion, for a maximum shaking weight of 7.5 kg.
- Electronic adjustment of speed and timer
 - LED display for speed and time adjustment
 - Wide range of attachment combinations makes it possible to use almost all shapes and sizes of vessels
 - Attachments are not included in delivery, please order separately

| Technical Data | | | |
|--|--------------------------------|-----------------------------------|-----------------------|
| Type of movement [-] | reciprocating | Orbital diameter [mm] | 20 |
| Permissible shaking weight (incl. attachment) [kg] | 7.5 | Motor rating input [W] | 45 |
| Motor rating output [W] | 10 | Permissible ON time [%] | 100 |
| Speed min (adjustable) [1/min] | 20 | Speed range [1/min] | 0 - 300 |
| Speed display [-] | LED line | Timer [-] | yes |
| Timer display [-] | LED line | Time setting range [min] | 5 - 50 |
| Operating mode [-] | timer and continuous operation | Shaking table lock [-] | no |
| Dimensions (W x H x D) [mm] | 360 x 100 x 420 | Weight [kg] | 8.5 |
| Permissible ambient temperature [°C] | 5 - 50 | Permissible relative moisture [%] | 80 |
| Protection class according to DIN EN 60529 [-] | IP 21 | RS 232 interface [-] | no |
| Analog output [-] | no | Voltage [V] | 220 - 240 / 100 - 120 |
| Frequency [Hz] | 50/60 | Power input [W] | 45 |